



## Ultrasonic Skin Scrubber ANLAN ALCPJ07-02 (white)

Ref: 6953156299221

Ultrasonic Skin Scrubber ANLAN ALCPJ07-02 (white)

### Cavitation peeling ANLAN ALCPJ07-02

Do you dream of beautiful, smooth skin? Cavitation peeling by ANLAN will help you get the desired results! It offers up to 4 modes of operation - positive ion export, negative ion import, EMS microcurrent massage and ultrasonic vibration. With its help you will effortlessly reduce the appearance of wrinkles, smooth your skin and make it extra firm. The device will also help you fight acne and blackheads, give your complexion a healthy, beautiful glow and improve circulation.

### Reliable cleansing

Get rid of blackheads and impurities. ANLAN's cavitation peeling offers ultrasonic vibrations with a frequency of up to 24,000 vibrations per minute. It also uses the action of positive ions. As a result, it enables effective, deep cleansing. It effectively eliminates excess sebum, makeup residues, dirt and blackheads. It also exfoliates dead skin. As a result, the skin becomes smooth, shiny and radiant!

### Discover the power of negative ions

Increase the effectiveness of your cosmetics. ANLAN cavitation peeling can also generate negative ions, which make your skin absorb nutrients much better. As a result, your complexion will be perfectly moisturized and get a healthy glow. Enjoy optimal refreshment and be sure to provide your face with the best!

#### EMS massage

Provide yourself with a relaxing massage and enjoy a wonderful lifting effect. EMS microcurrents will strengthen facial muscles and increase skin elasticity, visibly reducing the appearance of wrinkles. They will also improve circulation and collagen production. As a result, your complexion will become smooth, firm and radiant.

#### Cordless design

Say goodbye to tangled cables and get more freedom of movement when performing skin care treatments. The ANLAN skin scrubber works wirelessly - it's equipped with a durable 500mAh battery. It stands out for its long working time, and it only takes about 2 hours to charge. You will also find a USB cable in the package, with which you can connect the device to the charger.

#### Easy to use

ANLAN's cavitation peeling will work well in many applications - you can successfully use it for skin care on the forehead, around the eyes, on the nose and on the cheeks. Its operation is very simple. First wash your face, then use a sauna or hot water to open your pores. Use the ultrasonic vibration mode to cleanse your skin, then apply your favorite skin care products and turn on the positive/negative ion mode. Finally, choose an EMS massage.

#### Thoughtful design

Skin care with ANLAN cavitation peeling is extremely comfortable. The device fits well in the hand and is easy to use. Its spatula is made of high-grade 304 steel and is curved at a 45° angle, so it glides smoothly over the skin and does not cause irritation. Possible transportation of the device is facilitated by the included pouch. No longer do you have to give up facial care when you go on vacation!

#### In the box

Cavitation peeling

Storage pouch

Charging cable

User manual

Brand

ANLAN

Name

Ultrasonic Skin Scrubber

Model

ALCPJ07-02

Color

White

Rated voltage

3.7-5V

Charging time

Approx. 2h

Battery capacity

500mAh

Working time in standby mode

Approx. 50-60min

Automatic shutdown

After 5 minutes / mode

Material

ABS

Dimensions

Approx. 171x40x14mm

Weight  
Approx. 75g

Preço:

**€ 25.50**

Saúde e beleza, Skincare equipment