



## IPL Crystal Freezing Point Hair Removal InFace ZH-18E (whit

Ref: 6971308400660

IPL Crystal Freezing Point Hair Removal InFace ZH-18E (whit

IPL laser epilator with cooling InFace ZH-18E (white).

Get rid of hair and enjoy smooth skin for longer! IPL laser depilator with cooling InFace ZH-18E allows you to get rid of unwanted hair from your body. It uses cooling technology and powerful IPL pulsed light. It allows you to choose from 5 available intensity levels. It also has an intelligent skin detection sensor. In addition, it is equipped with a practical LCD display.

Perfectly smooth skin in 4 weeks

After just 4 weeks you can enjoy touchable and silky skin! The energy of the device is 10 J. Use the epilator regularly for about 3 weeks to see the first results.

Cooling technology

The product uses an innovative cooling technology that provides even better skin protection. The system effectively reduces the amount of heat released and protects the epidermis, making the removal of unwanted hair painless, without burning or redness.

Precise depilation of the entire body

Model ZH-18E is equipped with a special sensor that intelligently detects the skin. If the device detects that the head is in contact with the skin, it will automatically perform a quick flash, without the need to press a dedicated button. Thanks to this solution, you can easily

remove unwanted hair even from hard-to-reach areas.

#### Tailored to different users

The epilator has 5 intensity levels, so you can perfectly match the intensity of the light to your skin tone or hair color. Levels 1 / 2 / 3 will work well for those with dark skin and black hair. While levels 4 / 5 are suitable for light complexions and light brown or black hair. The device is least effective for very light, red or gray hair - they contain too little melanin to absorb the light emitted in front of the epilator.

#### Instant hair removal

It only takes about 10 minutes to get rid of hair from the whole body! The InFace ZH-18E laser epilator works as much as 2 times faster than standard laser epilators - one flash takes just 0.8s. What's more, the device emits light with a wavelength in the range of 510 - 1200 nm, which penetrates hair roots and weakens hair follicles.

#### Easy to use

Operating the device will not cause you any problems. To begin, epilate the area of your body where you want to use the device. Then connect the product to the adapter and put on the protective goggles you will find in the kit. Set the appropriate values and bring the epilator close to your skin to start removing hair. We recommend starting epilation with level one, and gradually increasing the level on subsequent treatments.

#### Included

InFace laser epilator

AC adapter

razor x3

protective goggles

instruction manual

Manufacturer

InFace

Model

ZH-18E

Color

White

Light

IPL

Exposure levels

5

Wavelength

510~1200 nm

Cooling

10°C

Fan speed

6600 RPM

Irradiation energy

10 J

Operating current

4A

Operating voltage

12 V

Preço:

Antes: € 150.9948

**Agora: € 126.00**

Saúde e beleza, Equipamentos para cuidado corporal