



Phototherapy Hair Removal Light ANLAN ALTMY02-US02 with US p

Ref: 6953156299160

Phototherapy Hair Removal Light ANLAN ALTMY02-US02 with US p

Laser epilator ANLAN ALTMY02-US02 with US connector

Say goodbye to old hair removal methods! With the Anlan device you will get rid of them in just 8 weeks. After just 5 weeks the hair becomes visibly thinner, and after 8 weeks you will notice hair loss and no new hair formation. The US02 is equipped with two safe sensors - temperature and skin. In addition, it has 5 levels of energy settings. Adjust the level to suit your body part.

Do away with traditional epilation methods

You no longer need to spend masses of money on beauty salon visits. Now you can quickly and painlessly get rid of your hair without leaving home! It only takes 4 minutes to enjoy smooth legs in 8 weeks. For the lip or bikini area, it only takes 1 minute. What's more, you can easily charge it using the charger with US plug and adapter included.

Enjoy smooth skin in 8 weeks

The device reaches a peak power of 5.3 J / cm². After 1-4 weeks, hair growth becomes gradual. By week 5-8, on the other hand, the hair starts to become thin and soft. After just 8 weeks, hair starts to fall out and no new hair appears. Note: When using the device for the first time, it is recommended to start at level 1 and gradually increase the power with use.

Equipped with two safety sensors

The temperature of the head constantly increases with the radiation, and the temperature sensor the unit is equipped with is designed to

block radiation until the main head cools down. To prevent erroneous radiation, it has a skin sensor that detects its adhesion. If the light indicators are not in full contact with the skin, they will not light up.

5 levels of energy settings

Adjust the energy level to the body part of your choice. If you are using the device for the first time, start at level 1 and increase with use. The highest level gives the best results. This will allow you to adapt the operation of the device to your needs.

Convenience of use

Before using the device, remove all hair from the areas you want to expose. Be sure to wear glasses to avoid harmful radiation. Adjust the energy level to your needs, it is recommended to start with the lowest. Hold the product vertically to the epilation area and then press the flash button. Afterwards, don't forget to moisturise your skin.

Included

Depilator

Adapter US — EU

Shaver

Glasses

User manual

Case

Specification:

□ Producer

□ Anlan

□ Model

□ ALTMY02-US02

□ Wavelength

□ 530-1150 nm

□ Energy level

□ 1.3-5.3 J/cm²

□ Optional energy level

□ 5

□ Treatment area

□ 4 cm²

□ Power supply

□ 24 V / 2 A

□ Speed

□ 1.4-3 seconds/pulse

□ Maximum light pulse

□ 500000

```
.table-wrapper{width:100%;overflow-x:auto}.table-wrapper>table{border-spacing:0}.table-wrapper>table>tbody>tr>td,.table-wrapper>table>tbody>tr>th{padding:10px 35px;text-align:left}.table-wrapper>table>tbody>tr>th{vertical-align:top}.table-wrapper>table>tbody>tr:nth-child(even){background-color:#f7f7f7}
```

Preço:

€ 42.50

Saúde e beleza, Equipamentos para cuidado corporal