



## Electric EMS neck massager ANLAN ALMJY01-02

Ref: 6953156299863

Electric EMS neck massager ANLAN ALMJY01-02

Electric neck massager by ANLAN

Enjoy smooth, firm and visibly rejuvenated skin. The neck massager by ANLAN offers several modes for comprehensive skin care. It allows deep cleansing and can provide EMS massage. It affects the skin with 45°C temperature and light in 3 colors - red, blue and green. It also uses ion action. It is equipped with a durable 500 mAh battery.

3 working modes

The 3 available modes of the massager can meet the diverse needs of users. The EMS mode allows electrostimulation of muscles and helps improve collagen production. It smooths wrinkles and increases skin elasticity. Clean mode cleans pores and adds radiance to the complexion. Hot mode uses a temperature of 45°C to gently stimulate and relax the skin, so it contributes to skin rejuvenation.

Light therapy in 3 colors

The massager also offers light therapy in 3 colors. Red light removes wrinkles, regenerates damaged skin and improves collagen production. Blue light gets rid of acne and soothes inflammation. Green light promotes blood microcirculation in the skin and reduces swelling. All this will allow you to improve the appearance and condition of your complexion without having to spend money on visits to a beauty salon.

#### Sonic vibration and 45°C temperature

With the help of a massager from ANLAN, you can also improve circulation and slow down skin aging. The device offers a special massage using sonic vibrations and a high temperature of 45°C. In this way, it smooths the skin, restores its elasticity and eliminates swelling.

#### Take advantage of the effect of ions

Positive and negative ions allow you to effectively cleanse your pores, and make your skin better absorb nutrients. This allows you to get rid of acne and solve the problem of dry, dull skin. See how little it takes to make your skin smooth, firm and radiant!

#### Durable rechargeable battery

The cordless massager lets you free yourself from wires and enjoy almost unlimited freedom of movement. Its long-lasting 500 mAh battery ensures a long run time. Charging the device takes only about 2 hours. To connect the massager to a power source, use the included USB cable.

#### Multiple uses

This practical massager will prove useful in many applications. With its help you can nurture not only your neck, but also your neck, forehead and cheeks. This will comprehensively care for the appearance and condition of your skin. Give yourself a relaxing massage, deep cleansing and effective light therapy to delight in a smooth, firm complexion.

#### Included

Neck massager

USB charging cable

User manual

□Manufacturer□

□ANLAN

□Name□

□Neck Beauty Device

□Model□

□ALMJY01-02

□Material□

□ABS

□Output voltage□

□3,7 V

□Battery capacity□

□500 mAh

□Charging time□

□Approx. 2 h

□Power□

□2,5 W

□Input voltage□

□5 V 1 A

□Dimensions□

□89 x 47 x 120 mm

Preço:

**€ 34.50**

Saúde e beleza, Massage Guns