











Ref: 6953156299276

Intelligent eye massager ANLAN ALMYY06-02 Intelligent eye massager ANLAN ALMYY06-02

Intelligent eye massager ANLAN ALMYY06-02

Give yourself a relaxing massage at home. Regular use will make your skin relaxed and smooth. The device has 3 types of light - red, yellow and blue. Are you struggling with lack of elasticity and sagging skin? With their help, you can repair and make your skin more elastic. In addition, thanks to the adjustable temperature (38 - 45 $^{\circ}$ C) you can eliminate dark circles under your eyes. The massager is equipped with a slanted 45 $^{\circ}$ design, for better comfort of use.

Thermotherapy

The high vibration frequency of 11500 2000 rpm provides a unique experience. The device smooths out fine lines and swelling in just 20 days. Regular use stimulates circulation and nourishes, thus lifting the complexion, brightening the eye area and giving a feeling of relaxation.

Overcome your problems once and for all

The massager helps nutrients to be absorbed much deeper and better than we would do manually. You can say goodbye to dark eye circles, wrinkles or bags under your eyes.

Tri-color light therapy

Mobilemax - Gadgets and Accessories Rua Vasco da Gama n61, 8005-411, Estoi - Faro, Portugal (Continente) Telefone: +351962484153 (Chamada para rede móvel nacional) Email: info@mobilemax.pt



Using red light you will improve cell function, collagen stimulation will increase, thus skin will become more elastic and soft. Blue light kills bacteria, it works well in treating acne. You will get a soft and shiny complexion thanks to proper distribution of melanin, which yellow light will help you with.

Easy to use

Before using the massager, cleanse your face thoroughly and select the mode you are interested in. Move the roller back and forth making small circles from the center to the lower eyelid. Move from the inside of the lower eye to the inside of the upper eye. You can use the product on your entire face, including your lips.

Use wherever you need

The device can be in standby mode for up to 3 hours. However, it will automatically turn off after 15 minutes. The angled design at 45 allows for a perfect fit to the skin and better massage results.

Included

Massager

Charging cable

Case

User manual

Brand

Anlan

Model

ALMYY06-02

Standby time

3 hours

Auto power off

15 minutes

Charging time

2 hours

Vibration frequency

11500 2000 rpm

Temperature

38 - 45 °C

Size

130 x 22 x 22 mm

Weight

42 g

Battery Capacity

400mAh

Power

2.5 W

Power supply

5V 1A

Nominal voltage

3.7V

Preço:

€ 29.00



Telefone: +351962484153

Email: info@mobilemax.pt

(Chamada para rede móvel nacional)

Saúde e beleza, Massage Guns